## Stimulator Accessories & Services



CS Medical Systems offers a full line of accessories and replacement parts for our Interferential and Galvanic Stimulators. Our accessories will protect your investment and keep your stimulators running in optimal condition.

Call customer service or your representative for more information.

#### **Carrying Pouch/Protector**

The Carrying Pouch/Protector is a padded neoprene pouch designed to protect the stimulator from scratches. The pouch features a hook & loop closure, open access to the lead

terminals, and a convenient belt clip that allows the patient to attach the stimulator to their belt while in use. Order part number 00617A. For CS6101 Stimulators only.



# Repair & Rechargeable Battery Replacement Services

CS Medical Systems offers full repair services for our line of Interferential and Galvanic Stimulators. In order to provide maximum performance, the CS6101 Stimulator is calibrated to its rechargeable battery at our factory. For this reason,

when replacing the battery pack, the unit must be sent to Control Solutions. The unit will be fully tested and calibrated when the battery is replaced.



### **Lead and Electrode Sets**

Replace old and worn lead and electrode sets with our direct factory replacement lead sets.

Lead Order Number	Use with Stimulator
CS2110B	CS6101, CS6103
CS2112	CS6102
Electrode Order Number	Use with Stimulator
CS1304	CS6101, CS6103
CS1304A	CS6102





### **AC Adapters**

Our AC adapters save money in battery replacement costs. Two models are available for the different Stimulator models.

Order Number	Use with Stimulator
CS1210A	CS6101, CS6103
CS1211	CS6102









BR 0043-00 (012112)





CS Medical Systems c/o Control Solutions LLC 2520 Diehl Road, Aurora, IL 60502 Tel: 630.806.7062 Fax: 630.806.7065 www.csmedsys.com